

Long-Term Modeling and Monitoring of Neuromusculoskeletal System Performance Using a Novel Skin-Compliant Epidermal Sensing System

Kaiwen Yang

Skin-Compliant Epidermal Sensing System (ESS)




Figure 1: Tattoo-Like Sensor Geometry




Figure 2: Cut-and-Print Manufacturing Process of New Sensor





Figure 1: Experiment Set Up


Experiment Types:

- **Same-day repetitive fatigue trial**
2 mins for per trial X 8 ,30 mins rest time, one person
- **Multi-day repetitive fatigue trial**
2 mins per trial, 2 trials per day in 3 consecutive days, one person
- **Repetitive Fatigue and recovery trial**
Gripping trial :2 mins X 3, 90 mins rest time
Recover detection trial : 10 secs X 9, 10 mins rest time

System Input: EMG 1000HZ




System Output: Grip Force 100HZ




Trial Code: C2016_04_24_Luke

Data Processing on EMG signal



Raw EMG Signal


ARMAX Formulation



System Output : $F(t)$


System Input : $\vec{X}(t) = \begin{bmatrix} f^0|t \\ f^1|t \\ f^2|t \\ E \end{bmatrix}$

$$F(t) = \sum_{n=1,2,\dots,N} \phi_n F(t-n) + \sum_{n=1,\dots,N-1} \theta_n a(t-n) + \vec{X}(t)$$



Modeling window: 20 seconds


WHY??



Model Order (11,12) for trial code:
'C2016_04_24_Luke'


- 1) Set Initial ARMAX order $(N, N-1) \leftarrow (2, 1)$
- 2) **for** $i \in \{1, 2, 3\}$ **do**
- 3) **if** ARMAX($2 \times i, 2 \times i - 1$) vs ARMAX($2 \times i + 2, 2 \times i + 1$) can fix model order **then**
- 4) Model order fixed \leftarrow **True**
- 5) **else**
- 6) Model order fixed \leftarrow **False**
- 7) **end if**
- 8) **end for**
- 9) **if** \sim Model order fixed **then**
- 10) Compute AIC of model up to order $N=15$
- 11) Find \tilde{n} s. t. : $\operatorname{argmax}\{AIC(\tilde{n}) - AIC(\tilde{n} - 1)\}$
- 12) Fix ARMAX model order as ARMAX $(\tilde{n}, \tilde{n} - 1)$
- 13) **end if**
- 14) Check autocorrelation plot for adequacy

Model Evaluation



Lag vs Residual Correlation


Model is OK consider a system sampled at 1000 HZ with extremely noise EMG input.



Poles Location

$-0.1708 + 0.8983i$
 $-0.1708 - 0.8983i$
 $-0.7420 + 0.0000i$
 $-0.6586 + 0.0000i$
 $0.4710 + 0.6497i$
 $0.4710 - 0.6497i$
 $1.0000 + 0.0000i$
 $0.9461 + 0.0409i$
 $0.9461 - 0.0409i$
 $0.4242 + 0.0000i$
 $-0.0407 + 0.0000i$


Fatigue Tracking





$$Fatigue\ Index\ (FI) := D_{KL} (P||Q_T) = \sum_i P(i) \cdot \ln \frac{P(i)}{Q_T(i)}$$

FI measures the similarity between current system condition and initial reference system.
 FI close to zero means current system is similar to unfatigued system.
 FI close to one means current system is different with unfatigued system.

Results







Conclusion:

High nonlinearity and noisy input and output causes residual correlation not close to zero.

By using an ARMAX Model we are able to model unfatigued human forearm neuromuscular dynamics and monitor fatigue as gripping continue.

Future Work:

Nonlinear neuromusculoskeletal system modeling with joint dynamics.